# HOFFMANN'S CARBOXYLAT CEMENT

Zink-Polycarboxylat Cement

---

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Zinc polycarboxylate cement

---

## Instruction for Use for the dentist

### Composition:

powder contains: zinc oxide, magnesium oxide;

liquid contains: polyacrylic acid

### Areas of Application:

Fixation of crowns and inlays, cavity linings

### Dispensing Instructions:

Cavity lining consistency: 3.6 g powder : 1.0 g liquid

Fixation consistency: 2.9 g powder : 1.0 g liquid

Dispensing and testing climate for testing in accordance with DIN EN ISO 9917:

Mix quantities indicated at  $23 \pm 1^\circ\text{C}$  and a relative humidity of  $50 \pm 10\%$ .

The dispensing instructions also apply when mixing the liquid with HOFFMANN'S UNIVERSAL POWDER.

### Processing Recommendations:

A clean, dry glass plate at ca.  $20^\circ\text{C}$  ( $68^\circ\text{F}$ ) and a non-discolouring spatula are recommended for mixing.

### Mixing:

Divide the measured quantity of powder into 4 portions (1/2 1/4 1/8 1/8). Beginning with the largest quantity, stroke the powder into the measured liquid within 60 seconds. The mixture should appear thicker as usual for zinc phosphate cement, yet it is just as easy to work. The filling consistency has been reached when on lifting the spatula, the peak drawn out forms a small hook and does not sink back into the paste.

The fixation consistency has been reached when the peak formed on lifting the spatula slowly falls back into the paste.

Discard remaining powder or impure liquids.

**Mixing time:** 60 seconds

**Processing time (at  $23^\circ\text{C}$ ,  $72^\circ\text{F}$ ):**

For fixation: 2 minutes

For cavity linings: 1 minute

**Setting Time from end of mixing:** 2.0 - 4.0 minutes.

### Special Information:

- The mixed cement should never be applied directly onto dentine near to pulp or onto exposed pulp.
  - Apply the thinnest possible layer for fixations.
  - The cement must be applied to the dry tooth surface and must set without contact to moisture.
  - Remove excess material from neighbouring teeth, facets etc. immediately.
  - Do not expose liquid to sunlight and frost.
  - Close containers immediately after removing material.
  - Clean instruments and glass plate with water before the cement sets.
- Zinc polycarboxylate cement can be dissolved using a cement-solvent after setting.

### Colours:

1 (white), 2 (bluish white), 3 (yellowish white), 4 (light yellow), 5 (yellow), 6 (golden), 7 (golden brown), 8 (pearl grey), 9 (grey), 10 (greenish grey), 11 (bluish grey), 12 (brown), 13 (greyish brown), 14 (light pink), 15 (pink)

### Package Sizes:

Powder: 100 g, Liquid: 40 ml,

Set: 35 g powder and 15 ml liquid

### HOFFMANN'S CARBOXYLATE CEMENT

- Do not use after the expiration date
- Compliant with DIN EN ISO 9917.